

9. Aldol Reaction

M. Jones: Condensation Reactions, Aldol reaction, Chapter 17, Section 3, pgs 840-850.

This procedure has been adapted from the microscale procedure described in the third edition of *Macroscale and Microscale Organic Experiments* by Kenneth L. Williamson (Houghton Mifflin, Boston, 1999).

Background

Aldol Reaction. In the experiment, you will perform a base-catalyzed, condensation reaction using benzaldehyde and acetone (see Figure 1). A condensation reaction is one, which condenses two or more molecules to make one single compound.




Figure 1. The overall reaction.

First let's define some different types of compounds. An aldehyde is a molecule with an R-HC=O structure. There is a carbonyl group (C=O). The carbon of the carbonyl group is flanked by a hydrogen (H) and an R group (R does not equal H). A ketone is a molecule with an R₂C=O. The R groups do not have to be the same. In this reaction, two molecules of benzaldehyde (aldehyde) are condensed with one molecule of acetone (ketone). Acetone acts as a nucleophile which adds to the carbonyl carbon of benzaldehyde. The **first step** in the reaction is to form the nucleophile, which is an enolate ion (Figure 2). This is a deprotonation step where a hydroxide ion pulls off a proton from the alpha carbon (carbon adjacent to the carbonyl carbon) to yield a resonance-stabilized enolate. The nucleophile adds to the carbonyl carbon of benzaldehyde in **step 2**. The resulting alkoxide

ion is protonated in **step 3** to form the "true" Aldol product which has both alcohol (OH) and carbonyl (C=O) functionalities.




Figure 2. General reaction mechanism for the condensation of one molecule of benzaldehyde with one molecule of acetone.

With heating, this product eliminates water (dehydration) to form an α,β -unsaturated ketone. This happens first by a deprotonation step (**step 4**) with sodium hydroxide to form a resonance-stabilized carbanion. Then in **step 5**, a hydroxide ion is eliminated to form the α,β -unsaturated ketone called benzalacetone (mp 42 °C). The entire reaction sequence is repeated to condense another molecule of benzaldehyde to the second alpha carbon of acetone and form dibenzalacetone.

Cautions and tips:

- Make sure all glassware is clean.
- Before you start the experiment, start the water bath and put ethanol on ice to save time.
- If no crystals form, add a few drops of water to the solution and scratch the side of the tube with a stir rod.

Experiment

Prepare an ice water bath for the recrystallization. Combine 0.6 mL of benzaldehyde (PhCOH), 0.2 mL of acetone (CH₃COCH₃), 8 mL of ethanol (CH₃CH₂OH), and lastly 4.0 mL of 5 M aqueous sodium hydroxide solution into your reaction tube and cap the tube. Make sure to record all amounts of starting materials. Shake the tube a few times every minute for a total of 30 minutes. After that time period, remove the liquid from the crystals using a pipette. Add 10 mL of water and shake the tube. Remove the water with a pipette and wash the tube two more times with water. Collect the crystals by vacuum filtration. Recrystallize the resulting crystals using warm ethanol. Vacuum filter and wash the crystals with a cold solution of ethanol/water (70:30). Dry and weigh the crystals, take a melting point and report the percent yield.